Technovit® 4000

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GB User instructions Technovit 4000

Application

For metallographic analyses:
gap-free embedding of ground sections

Material

Quick-hardening, cold-curing three-component resin based on modified polyester in the form of powder and syrup.

Properties

After mixing, Technovit 4000 is pourable for approximately 4 min; its good flow properties ensure that even the finest details are filled. Curing time at room temperature (approx. 22 °C): 6 – 13 min.
When correctly processed, the material's very good adhesion to metal and extraordinarily low shrinkage rate of 0.1–0.2 % allow gap-free embedding of even complicated specimens. Good abrasion resistance combined with excellent marginal seal guarantee precision-cast edges.
After polymerisation, Technovit 4000 can be mechanically processed (for example, ground, cut, drilled, polished).

Solubility

Technovit is resistant to acids and alkalis; when cured, it is insoluble in esters, ketones and chlorinated hydrocarbons. Polyethylene and Hostaphan films are not subject to chemical attack.

Processing

Preparations

The surfaces of the specimen which come into contact with Technovit 4000 must be clean and degreased with acetone. Polyethylene moulds should be used for embedding specimens; use of the Kulzer embedding moulds is recommended.

Mixing

Mix two parts by weight or two parts by volume of syrup I with one part by weight or one part by volume of syrup II. Approx. two to three parts by weight or by volume of powder are then stirred into the syrup mixture.

Examples of dosing

www.kulzer-technik.com → Metallography → cold embedding

By weight or by volume

two parts of syrup I
one part of syrup II
two parts of powder

or

two parts of syrup I
one part of syrup II
three parts of powder

Syrups I and II can be pre-mixed to facilitate processing, but should not be stored for more than 5 days. In such cases, processing is as with a two-component material. (The powder is stirred into the pre-mixed syrups in the ratio 3:2 to 1:1 liquid/powder).

The powder must be thoroughly stirred in to ensure homogeneous, lump-free distribution.

Thereafter, let the dough soak for approx. ½ minute.

Avoid a whipping motion when mixing, as this would lead to the formation of air bubbles.

At least mix dough carefully and slowly pour a thin stream of the mixture over the specimen, to allow any air in the material to escape. Mixing vessels and spatulas can easily be cleaned with esters, ketones and chlorinated hydrocarbons as long as the resin is still plastic; once hardened, Technovit 4000 is very difficult to remove mechanically.

Curing

Technovit begins to cure when the powder is added. When mixing larger quantities, the polymerisation heat can be carried off by placing the vessel in cold water (do not immerse).

Devesting

Because of the good adhesive properties and the low level of contraction, the deformation is recommended after the polymerisation while the material is still warm.

Note

Technovit 4000 should be stored in a cool, dark place, in tightly sealed containers. Contamination must be avoided.

Delivery units

Combination packs

Technovit 4000

750 g powder
500 ml sirup I
250 ml sirup II
measuring spoon

1500 g powder
1000 ml sirup I
500 ml sirup II
measuring spoon

Single package

Technovit 4000

1500 g powder
1000 ml sirup I
500 ml sirup II

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.
www.kulzer-technik.com

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

* For more information please visit www.kulzer-technik.com

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® = marchio registrato, Kulzer GmbH, Hanau
KULZER GmbH